

# Y-maze RBRCLA\_YMZ\_001

## Purpose

The Y-maze is used to assess spontaneous alternation performance as an index of active retrograde working memory in rodents.

## Experimental Design

- **Minimum number of animals** : 7M + 7F
- **Age at test**: Week 51

## Equipment

Y shaped arena with three arms of identical length at 120° angles.

Camera and computer with tracking software if using video tracking.

## Procedure

1. Transport animals to testing room at least 15 minutes prior to testing.
2. Ensure lighting conditions are as desired and if using video tracking that the camera is working.
3. Remove mouse from its cage and place in the start arm. Allow the mouse to explore the arena for 8 minutes.
4. Either manually score the entries or use the automatic tracking while the mouse explores the arena for the duration of the test. Arm entry is defined differently depending on the scoring system.
  - a. Manual: arm entry is scored when all four paws enter the new location.
  - b. Automated tracking: arm entry is scored when the centre of the mouse enters the new location.
5. Remove mouse from the arena and place back in home cage.
6. Clean area with consistent disinfectant (i.e. ethanol/clidox) before testing next mouse.

## Notes

Minimum of 5 arm entries are required for the ratios to be calculated to ensure the accuracy of the derived values.

## Parameters and Metadata

## Tracking equipment manufacturer RBRCLA\_YMZ\_025\_001 | v1.1

procedureMetadata

Req. Analysis: true

Req. Upload: false

Is Annotated: false

Options: Noldus, Accuscan, O'hara Co.Ltd,

---

## Total arm entries RBRCLA\_YMZ\_002\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

---

## Same arm entry ratio per total arm entries RBRCLA\_YMZ\_009\_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: %

Derivation:

ifElse(greaterThan('RBRCLA\_YMZ\_002\_001',4), mul(div('RBRCLA\_YMZ\_008\_001', 'RBRCLA\_YMZ\_002\_001'), 100),")

---

## Tracking method RBRCLA\_YMZ\_017\_001 | v1.1

procedureMetadata

**Req. Analysis:** true

**Req. Upload:** true

**Is Annotated:** false

**Options:** None, Video,

-----

**Maze model** RBRCLA\_YMZ\_023\_001 | v1.0

procedureMetadata

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** false

**Options:** JAX\_YM, YM-3002,

-----

**Tracking equipment model** RBRCLA\_YMZ\_026\_001 | v1.0

procedureMetadata

**Req. Analysis:** false

**Req. Upload:** false

**Is Annotated:** false

**Options:** Opto-Varimex 4, Time YM2,

-----

**Latency to leave start arm** RBRCLA\_YMZ\_001\_001 | v1.0

simpleParameter

**Req. Analysis:** false

**Req. Upload:** true

**Is Annotated:** true

**Unit Measured:** s

---

## Tracking equipment ID RBRCLA\_YMZ\_024\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Options: Ethovision, Y tracker, YM-1020,

---

## Number of spontaneous alternations RBRCLA\_YMZ\_004\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

---

## Software RBRCLA\_YMZ\_028\_001 | v1.1

procedureMetadata

Req. Analysis: true

Req. Upload: false

Is Annotated: false

Options: Ethovision 13.0, Time YM2 software,

---

## Light level RBRCLA\_YMZ\_013\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: Lux

Options: 50, 20, 4,

---

## Number of triplets RBRCLA\_YMZ\_003\_001 | v1.2

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Derivation:

ifElse(greaterThan('RBRCLA\_YMZ\_002\_001', 1), sub('RBRCLA\_YMZ\_002\_001', 2),0)

---

## Duration at centre prior to arm re-entry RBRCLA\_YMZ\_019\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: s

Options: 1,

---

## Experimenter ID RBRCLA\_YMZ\_020\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

---

## Cleaning agent/disinfectant RBRCLA\_YMZ\_016\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Options: 70% Ethanol solution, 6% sodium hypochlorite solution,

---

## Alternate arm entry ratio per number of triplets RBRCLA\_YMZ\_007\_001 | v1.3

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: %

Derivation:

ifElse(greaterThan('RBRCLA\_YMZ\_002\_001',4), mul(div('RBRCLA\_YMZ\_006\_001', sub('RBRCLA\_YMZ\_002\_001', 2)), 100),")

---

## Number of alternate arm entries RBRCLA\_YMZ\_006\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Maze ID RBRCLA\_YMZ\_021\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Options: Arena 1, Arena 2, y-maze 1,

---

## Arm dimensions RBRCLA\_YMZ\_012\_001 | v1.1

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Unit Measured: cm

Options: 33.65 L x 6 W x 15 H, 40 L x 3 W x 12 H,

---

## Alternation ratio per number of triplets RBRCLA\_YMZ\_005\_001 | v1

.2

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: true

Unit Measured: %

**Derivation:**

ifElse(greaterThan('RBRCLA\_YMZ\_002\_001',4), mul(div('RBRCLA\_YMZ\_004\_001', sub('RBRCLA\_YMZ\_002\_001', 2)), 100),")

**Number of same arm entries by number of doublets** RBRCL

A\_YMZ\_008\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

**Alternating start arm** RBRCLA\_YMZ\_029\_001 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Options: Yes, No,

**Maze manufacturer** RBRCLA\_YMZ\_022\_001 | v1.0

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Options: The Jackson Laboratory, O'hara Co.Ltd,



---

## Acclimatisation time (minimum) RBRCLA\_YMZ\_011\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: min

Options: 60, 15, 30,

---

## Scoring method RBRCLA\_YMZ\_018\_001 | v1.1

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Options: Automated, Manual,

---

## Feacal boli (count) RBRCLA\_YMZ\_010\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Placement in start arm RBRCLA\_YMZ\_015\_001 | v1.1

procedureMetadata

Req. Analysis: true

Req. Upload: true

Is Annotated: false

Options: End of arm, Mid-arm,

-----

**Date tracking equipment last calibrated** RBRCLA\_YMZ\_027\_001 |  
v1.0  
procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

-----