Hot Plate HRWL_HOT_001

Purpose

To assess thermal sensitivity.

Experimental Design

Minimum number of mutant animals: 7 males + 7 females

Age at test: 9 weeks

Sexual dimorphism:

Procedure

- 1. Acclimatise the mice to the arena for 30 seconds on the day before testing.
- 2. On the day of the test, switch on the hot plate and bring up the temperature to 52° Celsius. Two people should be present for the test.
- 3. Set up the video recording camera to record a 360° view of the hot plate, using mirrors to reflect the images.
- 4. Place the mouse on the hot plate, depressing the foot pedal simultaneously.
- 5. Both experimenters observe the mouse from 2 different angles for any reaction to the stimulus.
 - If a reaction is observed then the foot pedal is depressed again and the mouse is removed from the hot plate.
 - The mouse is removed after 30 seconds, even if no reaction has been observed.
- 7. The latency to response and the type of response are manually recorded using the video output.

Notes

This procedure is a pilot study from the Pain Phenotyping Pilot

Parameters and Metadata

Body weight HRWL_HOT_001_001 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: true Is Annotated: false

Unit Measured: g			
Response latency simpleParameter	HRWL_HOT_002_001 v1.0		
Req. Analysis: false	Req. Upload: true	Is Annotated: false	
Unit Measured: s			
Type of response H simpleParameter	IRWL_HOT_003_001 v1.0		
Req. Analysis: false	Req. Upload: true	Is Annotated: false	
Options: Step, No reaction, Locomotion, Jump, Lick, Retropulsion, Shake, Foot roll,			
Equipment manufa procedureMetadata	cturer HRWL_HOT_004_	001 v1.0	
Req. Analysis: false	Req. Upload: true	Is Annotated: false	
Options: Bioseb,			

Equipment model HRWL_HOT_005_001 | v1.0

procedureMetadata

Req. Analysis: false	Req. Upload: true	Is Annotated: false		
Options: Dynamic Hot Plate,				
Maximum plate temperature HRWL_HOT_006_001 v1.0 procedureMetadata				
Req. Analysis: false	Req. Upload: true	Is Annotated: false		
Unit Measured: C				
Options: 52,				
Maximum time on hot plate HRWL_HOT_007_001 v1.0 procedureMetadata				
Req. Analysis: false	Req. Upload: true	Is Annotated: false		
Unit Measured: s				
Options: 30,				

Experimenter ID HRWL_HOT_008_001 | v1.0

procedureMetadata

Req. Analysis: false	Req. Upload: true	Is Annotated: false