## Electrocardiogram (ECG) UCDLA\_ECG\_003

### Purpose

To provide a high throughput method to obtain Electrocardiograms in a conscious mouse or an anesthetized mouse.

### **Experimental Design**

- Minimum number of animals : 5M + 5F
- Age at test: Week 55
- Sex: We would expect the results of this test to show sexual dimorphism

### Procedure

#### **Conscious Procedure with ECGenie**

- 1. The lead plates are to be snapped into place onto the top of the pre-amplifier tower. The covering is removed to reveal three gel coated pads surrounded by a sticking plate. The plate will need to be covered with the extra cover in the package.
- 2. Turn on the combined amplifier and the pre-amplifier tower.
- 3. Double click the icon ECG acquisition on the acquisition computer.
- 4. Open the ECG set up file (for default settings).
- 5. Place mouse on pad, lowering the Red Acrylic Cubby to surround the mouse on 3 sides discouraging escape.
- 6. Press Start.
- 7. After the desired acquisition time, (5-10 minutes) stop the reading. There will be one long reading.
- 8. Save the data.
- 9. For additional readings create a new session using the same settings as before.
- 10. When saving sections with good readings, highlight the selected area and then save.

#### **Anesthetized Procedure**

- 1. Place mouse in the induction chamber. Anesthetize with 2.5-4% isoflurane in oxygen.
- 2. Transfer mouse to a warmed platform and maintain with a nose cone at 2-2.5% isoflurane.
- 3. Fasten the mouse to the heated platform set to 37.5-39° Add electrode cream if using the foot surface electrodes on the platform surface. *Alternatively insert needle electrodes subcutaneously into the limbs using Lead Configuration I: Left Arm Right Arm (LA-RA) and/or Lead Configuration II: Left Leg Right Arm (LL-RA).*
- 4. Monitor the body temperature with a rectal probe if possible.

- 5. Open the recording software with appropriate settings and record an ECG for about 2 minutes.
- 6. Save the recording and allow mouse to recover.

### Notes

#### Data Analysis - Conscious Procedure with ECGenie

- 1. Open Emouse Analyses icon
- 2. Select ECG signals
- 3. Choose folder (all readings in folder will show)
- 4. Click PNN X (for mice: N-N> than <u>6</u> ms)
- 5. Choose file(s) by highlighting
- 6. Go
- 7. Bottom file is the corrected file
- 8. Red dots should be on peak of R waves, if image appears inverted click invert
- 9. Click Add, or minus if R waves are not marked with red dots or if too many are marked
  - L click to zoom in
  - R click to zoom out
- 11. "What if?" button to remove unwanted sections
  - L click image (zooms in)
  - L click left boundary
  - L click right boundary
- 13. Options- click more if want to exclude more sections
- 14. Undo available
- 15. Go
- 16. Here can input animal data if desired
- 17. Save- For the first mouse in in group, hit save, a new results folder will be created within the folder with the mouse data. Then can click quick save or next.
- 18. For the rest of the mice in the series, can hit quick save at this point- saves in last selected file will group all files together in same excel sheet.
- 19. Open Emouse Analyses icon
- 20. Select ECG signals
- 21. Choose folder (all readings in folder will show)
- 22. Click PNN X (for mice: N-N> than 6 ms)
- 23. Choose file(s) by highlighting
- 24. Go
- 25. Bottom file is the corrected file
- 26. Red dots should be on peak of R waves, if image appears inverted click invert
- 27. Click Add, or minus if R waves are not marked with red dots or if too many are marked
  - L click to zoom in
  - R click to zoom out
- 29. "What if?" button to remove unwanted sections
  - L click image (zooms in)
  - L click left boundary
  - L click right boundary
- 31. Options- click more if want to exclude more sections
- 32. Undo available
- 33. Go
- 34. Here can input animal data if desired

- 35. Save- For the first mouse in in group, hit save, a new results folder will be created within the folder with the mouse data. Then can click quick save or next
- 36. For the rest of the mice in the series, can hit quick save at this point- saves in last selected file will group all files together in same excel sheet



Examples of good readings

#### Data QC

Analysis room should be dim and quiet. Keep the door closed preferably while analysis is taking place.

Figure A. Taking a reading



Figure B. Saving a section of the reading



Figure C. Analysis phase, with the options to remove sections on the "What if?" button below.



#### Data Analysis - Anesthetized Procedure

- 1. Review the tracing and note any abnormal findings. Take images for upload where required.
- 2. Analyse all or a selection of the tracing (approximately 100 beats minimum).
- 3. If using LabChart software (ADInstruments), select Mouse preset and QTc = Bazett.
- 4. The software can automatically mark the P, QRS and T waves of each beat (Beats to average = 1) and calculate the interval and amplitude data. The average of all beats will be determined for each parameter.
- 5. Alternatively, an averaged tracing can be generated first (Block Averaging). Markers for the P, QRS and T waves are placed by the software but can be adjusted according to the centre criteria. Interval and amplitude data is then generated by the software.

Example of the LabChart8 ECG setting

Example of LabChart ECG Averaging View

#### LabChart ECG Averaging View



#### **Parameters and Metadata**

#### Number of signals UCDLA\_ECG\_001\_001 | v1.2

#### simpleParameter

Req. Analysis: falseReq. Upload: trueIs Annotated: false

\_\_\_\_\_

#### HR UCDLA\_ECG\_002\_001 | v1.1

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: true

Unit Measured: bpm

.....

<b>CV</b> UCDLA_ECG_003_00 <sup>-</sup> simpleParameter	1   v1.0	
Req. Analysis: false	Req. Upload: false	Is Annotated: true
Unit Measured: %		
<b>RR</b> UCDLA_ECG_004_00	1   v1.2	
simpleParameter		
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Unit Measured: ms		
<b>PQ</b> UCDLA_ECG_005_00 <sup>-</sup> simpleParameter	1   v1.0	
Req. Analysis: false	Req. Upload: false	Is Annotated: true
Unit Measured: ms		

#### **PR** UCDLA\_ECG\_006\_001 | v1.1

simpleParameter

Req. Analysis: false	Req. Upload: false	Is Annotated: true
Unit Measured: ms		
QRS UCDLA_ECG_007_( simpleParameter	001   v1.2	
Req. Analysis: false	Req. Upload: true	Is Annotated: true
Unit Measured: ms		
<b>ST</b> UCDLA_ECG_008_001 simpleParameter	v1.0	
Req. Analysis: false	Req. Upload: false	Is Annotated: true
Unit Measured: ms		

### **HRV** UCDLA\_ECG\_010\_001 | v1.0

simpleParameter

Req. Analysis: false	Req. Upload: false	Is Annotated: true
Unit Measured: bpm		
Mean SR amplitude	<b>Ə</b> UCDLA_ECG_012_001   v	1.1
Req. Analysis: false	Req. Upload: false	Is Annotated: false
Unit Measured: mV		
Mean R amplitude simpleParameter	UCDLA_ECG_013_001  v1. <sup>-</sup>	1
Req. Analysis: false	Req. Upload: false	Is Annotated: false
Unit Measured: mV		
rMSSD UCDLA_ECG_01 simpleParameter	4_001   v1.0	
Req. Analysis: false	Req. Upload: false	Is Annotated: true
Unit Measured: ms		

### **pNN5(6>ms)** UCDLA\_ECG\_015\_001 | v1.2

simpleParameter

Req. Analysis: false	Req. Upload: false	Is Annotated: false
Unit Measured: %		
Equipment ID UCDLA procedureMetadata	A_ECG_016_001   v1.0	
Req. Analysis: false	Req. Upload: true	Is Annotated: false

#### Equipment Manufacturer UCDLA\_ECG\_017\_001 | v1.0

procedureMetadata

Req. Analysis: true	Req. Upload: true	Is Annotated: false
---------------------	-------------------	---------------------

**Options:** Mouse Specifics, Inc., AD Instruments, World Precision Instruments, Indus Instruments,

### Equipment Model UCDLA\_ECG\_018\_001 | v1.0

procedureMetadata

Req. Analysis: true Req. Upload: true

Is Annotated: false

Options: ECGenie, ML870/p, ML826/FE132, Iso-DAM8A, ECGenie + gel pads, ML866, Mouse MonitorS, PowerLab: 4/35,

\_\_\_\_\_

#### Anesthetic UCDLA ECG 019 001 | v1.0

procedureMetadata

Req. Analysis: true	Req. Upload: true	Is Annotated: false

**Options:** Isoflurane, No anesthesia, Avertin, Tribromoethanol, Isoflurane: 1.5-2.5(%), Isoflurane: 1.5-3(%), Isoflurane: 2(%), Isoflurane: 2-2.5(%), Isoflurane: 1.5-2(%), Tribromoethanol: 250(mg/kg),

#### Experimenter ID UCDLA\_ECG\_020\_001 | v1.0

procedureMetadata

Req. Analysis: false	Req. Upload: true	Is Annotated: false

#### Noise level UCDLA\_ECG\_021\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

# Light level UCDLA\_ECG\_022\_001 | v1.0 procedureMetadata Req. Analysis: false Req. Upload: false Is Annotated: false \_\_\_\_\_ Date equipment last calibrated UCDLA\_ECG\_023\_001 | v1.1 procedureMetadata Req. Analysis: false Req. Upload: false Is Annotated: false \_\_\_\_\_ QTC UCDLA ECG 009 002 | v2.0 simpleParameter Req. Analysis: false Req. Upload: false Is Annotated: false Unit Measured: ms **Derivation:** div('UCDLA ECG 028 001',sqrt(div('UCDLA ECG 004 001',100)))

\_\_\_\_\_

#### Analysis Software UCDLA\_ECG\_024\_001 | v1.0

procedureMetadata

Req. Analysis: true Req. Upload: false Is Annotated: false

Options: eMouse, Matlab, PowerLab, LabChart,

.....

#### Waveform Image UCDLA\_ECG\_025\_001 | v1.0

seriesMediaParameter

Req. Analysis: false	Req. Upload: false	Is Annotated: false
Increments: Minimum 1		

#### Waveform Image Comment UCDLA\_ECG\_026\_001 | v1.0

simpleParameter

Req. Analysis: false	Req. Upload: false	Is Annotated: false
Description:		
Free text comment on the	waveform image. Use parame	eterAssociation of the image
parameter to link to this te	xt.	

#### Lead Configuration UCDLA\_ECG\_027\_001 | v1.0

procedureMetadata

Req. Analysis: false Re	<b>q. Upload:</b> false	Is Annotated: false
-------------------------	-------------------------	---------------------

Options: Lead I: LA-RA, Lead II: LL-RA,

#### **QT** UCDLA\_ECG\_028\_001 | v1.0

simpleParameter

Req. Analysis: false	Req. Upload: true	Is Annotated: true
Unit Measured: ms		
<b>QT Dispersion</b> UCDL simpleParameter	A_ECG_029_001   v1.0	
Req. Analysis: false	Req. Upload: false	Is Annotated: true
Unit Measured: ms		

### Data acquisition sampling rate UCDLA\_ECG\_030\_001 | v1.0

procedureMetadata

Req. Analysis: false	Req. Upload: false	Is Annotated: false
Unit Measured: kHz		
<b>Options:</b> 1, 2, 4, 10,		

### Abnormal ECG detected UCDLA\_ECG\_031\_001 | v1.0

simpleParameter

Req. Analysis: false	Req. Upload: false	Is Annotated: true
Options: Yes, No,		