OPT E9.5 IMPC_EOL_001

Purpose

To assess morphological defects in E9.5 embryos from lethal strains with OPT.

Experimental Design

- Set up timed matings with heterozygous mice
- Day 0 is defined as the midpoint of the prior dark cycle following the identification of a copulation plug.
- Collect embryos at E9.5 with 3 minimum homozygous and 1 littermate control
- Ensure embryos are bleed as much as possible.
- Reconstruct the scans according to Mark's instructions

Procedure

- 1. Set up timed mating with heterozygous animals. Dissect and collect >=3 homozygote embryos.
- 2. Process the embryos using the centre's staining and embedding protocol. Centre-specific methods are recorded within a document
- 3. Scan and generate reconstructions of embryos.

Notes

Embryos may be processed for Histopathology.

Parameters and Metadata

Embryo reconstruction IMPC_EOL_001_001 | v1.1

seriesMediaParameter

Req. Analysis: false Req. Upload: true Is Annotated: false

Increments: Minimum 1

Experimenter ID IMPC_EOL_002_001 | v1.0

procedureMetadata

Req. Analysis: false Req. Upload: true Is Annotated: false Equipment ID IMPC_EOL_003_001 | v1.0 procedureMetadata Req. Analysis: false Req. Upload: true Is Annotated: false Equipment Manufacturer IMPC_EOL_004_001 | v1.0 procedureMetadata Req. Analysis: false Req. Upload: true Is Annotated: false Equipment Model IMPC_EOL_005_001 | v1.0 procedureMetadata Req. Analysis: false Req. Upload: true Is Annotated: false

Image Pixel Size IMPC_EOL_006_001 | v1.0

Req. Analysis: false	Req. Upload: true	Is Annotated: false		
Unit Measured: μm				
E9.5 OPT preparation protocol IMPC_EOL_007_001 v1.0 procedureMetadata				
Req. Analysis: false	Req. Upload: true	Is Annotated: false		
Date equipment last calibrated IMPC_EOL_008_001 v1.0 procedureMetadata				
Req. Analysis: false	Req. Upload: false	Is Annotated: false		
Date of Sacrifice IMPC_EOL_009_001 v1.0 procedureMetadata				
Req. Analysis: false	Req. Upload: true	Is Annotated: false		

Req. Analysis: false	Req. Upload: true	Is Annotated: false		
Date of Scan IMPC_EOL_011_001 v1.0 procedureMetadata				
Req. Analysis: false	Req. Upload: true	Is Annotated: false		
Screenshots IMPC_EOL_012_001 v1.1 seriesMediaParameter				
Req. Analysis: false	Req. Upload: false	Is Annotated: false		
Increments: Minimum 1				