

# Body Weight RBRCLA\_BWT\_001

## Purpose

The body weight test measures the weight of the mouse in a time series, allowing monitoring of its evolution; also, it is required in many other procedures.

## Equipment

- Laboratory Balance

## Procedure

- Place a cage containing mice on an operation table
- Remove a mouse from its cage. Determine and record its weight.
- Return the mouse to the cage or continue with a procedure.

## Notes

The body weight of each mouse is to be measured at the 56th week. Optional additional weights may also be recorded.

## Parameters and Metadata

### Body weight RBRCLA\_BWT\_001\_001 | v1.0

simpleParameter

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Unit Measured: g

Description: body\_weight

-----

### General comments about the mouse RBRCLA\_BWT\_002\_001 | v1.0

simpleParameter

Req. Analysis: false                      Req. Upload: false                      Is Annotated: false

Description: general\_comments\_about\_the\_mouse

-----

**Equipment ID** RBRCLA\_BWT\_003\_001 | v1.0

procedureMetadata

Req. Analysis: false                      Req. Upload: true                      Is Annotated: false

Description: equipment\_name

-----

**Equipment manufacturer** RBRCLA\_BWT\_004\_001 | v1.0

procedureMetadata

Req. Analysis: false                      Req. Upload: true                      Is Annotated: false

Description: equipment\_manufacturer

Options: Ohaus, Denver Instrument, Scientech, Mettler Toledo, Kern & Sohn GmbH,  
Sartorius, Acculab, Radwag,

-----

**Experimenter ID** RBRCLA\_BWT\_005\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

Description: experimenter\_id

---

## Date equipment last calibrated RBRCLA\_BWT\_006\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: false

Is Annotated: false

---

## Equipment model RBRCLA\_BWT\_007\_001 | v1.0

procedureMetadata

Req. Analysis: false

Req. Upload: true

Is Annotated: false

**Options:** AV2101, SI-2002, SL-3100D, TP-202, PG3001-S, Adventurer Pro, 440-47N, AM100, XS802S, AB104-S, AV212C, PCB2000-1, 440-33, 440-47, 572-35, 770-14, EMB220-1, EMB500-1, EMB600-2, EW600-2M, BL310, Navigator 34120, AV213C, VIC-511, VIC-123, BP6100, CPA3202S, EMB 200-2, Entris 2201-1S, 201-10,

---

## Body weight curve RBRCLA\_BWT\_008\_001 | v1.1

seriesMediaParameter

Req. Analysis: false

Req. Upload: false

Is Annotated: false

Unit Measured: g

**Description:**

Derivation to collect all body weights measured (when available) and plot them correctly as time series.

**Derivation:** unimplemented("")

**Increments:** Minimum 1

-----