ECG (Electrocardiogram) (GMC) GMC_910_001

Parameters and Metadata

PR GMC_910_001_001 | v1.0

simpleParameter

Req. Analysis: false	Req. Upload: false	Is Annotated: true
Unit Measured: ms		
Description: PR		
Pdur GMC_910_001_002	v1.0	

simpleParameter

Req. Analysis: false	Req. Upload: false	Is Annotated: true
Unit Measured: ms		
Description: Pdur		

QRS GMC_910_001_003 | v1.0

simpleParameter

Req. Analysis: false

Unit Measured: ms

Description: QRS

QT GMC_910_001_004 | v1.0

simpleParameter

Req. Upload: false	Is Annotated: true
.0	
Req. Upload: false	Is Annotated: true
	.0

ST GMC_910_001_006 | v1.0

simpleParameter

Req. Analysis: false

Unit Measured: ms

Description: ST

RR GMC_910_001_007 | v1.0

simpleParameter

Req. Analysis: false	Req. Upload: false	Is Annotated: true
Unit Measured: ms		
Description: RR		

HR GMC_910_001_008 | v1.0

simpleParameter

Req. Analysis: false	Req. Upload: false	Is Annotated: true
Unit Measured: ms		
Description: HR		

Qamp GMC_910_001_009 | v1.0

simpleParameter

Unit Measured: mV

Description: Qamp

Ramp GMC_910_001_010 | v1.0

simpleParameter

Req. Analysis: false	Req. Upload: false	Is Annotated: true
Unit Measured: mV		
Description: Ramp		

Samp GMC_910_001_011 | v1.0

simpleParameter

Req. Analysis: false Req. Upload: false Is Annotated: true

Unit Measured: mV

Description: Samp

QTC GMC_910_001_012 | v1.0

simpleParameter

QRSamp GMC_910_001_013 | v1.0

simpleParameter

Req. Analysis: false	Req. Upload: false	Is Annotated: true
Unit Measured: mV		
Description: QRSamp		

Lead GMC_910_001_014 | v1.0

simpleParameter

Req. Analysis: false	Req. Upload: false	Is Annotated: false

Description: Lead
